

K20345C / K20745C / K21045C Freezers

USES: / CLINICAL APPLICATIONS / VACCINES & MEDICATIONS / BIOLOGICS / DIAGNOSTICS KITS / MOLECULAR BIOLOGY / INDUSTRIAL TESTING

K2's value performance laboratory / pharmaceutical chest freezers are designed for the most rigorous life science and commercial applications. These durable, easy to use units feature an intelligent microprocessor controller with digital display and adjustable temperature control.







K20345C K20745C K21145C

- Intelligent microprocessor controller with digital display and adjustable temperature control
- Visual and audible hi or low temperature, sensor failure, door ajar alarms
- Keyed door locks
- Magnetically sealed door with rapid cool-down mode after door opening for quick temperature recovery
- Cabinet constructed of white powder-coated 304 steel with CFC-free, high density insulation (2.4' thick)

- Easy swivel casters
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- One wire basket K20745C and K21145C
- Warranty: two-year full replacement for parts or compressor failure plus three additional years for compressor replacement

MODEL NUMBER	CAPACITY CU. FT. / liters	EXTERIOR Height x Width x Depth	TEMP RANGE (SETPOINT)	DOORS	SHELVES	DEFROST	CERTIFICATION	SHIPPING WT. LBS.
K20345C K20745C K21145C	3.5 / 100 7/198 11/311	32.5" x 27.5" x 20.5" 35.5" x 37" x 26.5" 35.5" x 50" x 26.5"	-10°C to -45C°C Preset -40°C	1 Lid	one basket	MANUAL	UL/cUL	102 lbs 113.5 lbs 176 lbs