

K205ULT / K210ULT / K221ULT / Ultra-low Temp Freezers

USES: / CLINICAL APPLICATIONS / VACCINES & MEDICATIONS / BIOLOGICS / DIAGNOSTICS KITS / MOLECULAR BIOLOGY / INDUSTRIAL TESTING

K2 ULT models feature space-saving, vacuum-sealed cabinet insulation designed to maximize space and cooling efficiency. K2 ULT's are equipped with intelligent, adjustable controls, alarms and data logger ports that make them ideal for frozen vaccine storage where temperature monitoring is critical. The upright ULTs feature adjustable stainless-steel shelves with both interior and exterior doors for maximum temperature stability.



K205ULT

K210ULT

K221ULT

K226ULT

- Intelligent microprocessor digital LCD display, password protected.
- Cabinet and ambient air temperature displays
- Adjustable temperature control
- Audible and visual high or low temperature alarms
- Power failure audible alarm
- Low battery audible alarm
- Sensor failure alarm
- Refrigerant: CFC-Free
- Shelves are adjustable in 1" inch increments
- Exterior and interior cabinet constructed of stainless steel with PUF-VIP, high density insulation
- Easy swivel casters
- Remote dry alarm contacts
- 0.5" temperature probe and data logger access port
- CE listed
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- Warranty: two-years full replacement for parts or compressor failure

MODEL NUMBER	CAPACITY CU. FT. / liters	EXTERIOR Height x Width x Depth	TEMP RANGE (SETPOINT)	DOORS	SHELVES	DEFROST	CERTIFICATION	SHIPPING WT. LBS.
K205ULT	5 / 140	39.25" x 32.25" x 28.75"	-25°C to -86°C		0			275 lbs
K210ULT	9.88 / 280	83.5" x 27.0" x 32.0"	-50°C to -86°C	1	3	AUTO	UL / cUL	305 lbs
K221ULT	20.84 / 590	75.60" x 36.22" x 30.71"	-50°C to -86°C		3			685 lbs
K226ULT	25.7 / 730	77.55" x 36.22" x 37.48"	-50°C to -86°C		3			725 lbs