

K212GDR Pharmaceutical / Laboratory

USES: / CLINICAL APPLICATIONS / VACCINES & MEDICATIONS / BIOLOGICS / DIAGNOSTICS KITS / MOLECULAR BIOLOGY / INDUSTRIAL TESTING

K2's 14 cu. ft. high-performance laboratory / pharmaceutical refrigerators are designed for the most rigorous life science, research, healthcare, clinical, and commercial applications. Featuring an intelligent microprocessor controller with digital display and adjustable temperature control, temperature stability of +2.0°C,,, and -2.0°C from the set point. Environmentally friendly, natural refrigerants are fully compliant with new EPA/ SNAP/UL regulations.



- Intelligent microprocessor controller with digital display adjustable by 0.1 degrees
- Temperature control preset +4°C
- Visual and audible hi or low temperature, senor failure, door ajar alarms
- Security keyed door locks
- Lighting
- Magnetically sealed door with rapid cool-down mode after door opening for quick temperature recovery
- Door-activated and exterior switch-activated interior LED lighting
- Cabinet constructed of white powder-coated 304 steel with CFC-free, high-density insulation (2.0' thick)
- R290 refrigerant, an environmentally green option compliant with future EPA and SNAP non-ozone depleting mandates

- Easy swivel casters
- Remote dry alarm contacts
- 0.75 inch temperature probe and, one data logger access port
- UL and cUL listed (61010-1)
- Energy Star high performance certified and listed
- Shelves are adjustable by .5 inch increments
- Built and shipped from ISO 9001:2015 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- Dual pane tempered glass door
- Warranty: two-year fullreplacement for parts or compressor failure plus three additional years for compressor replacement









| MODEL NUMBER | CAPACITY CU. FT. / liters | INTERIOR Width x Depth x Height | EXTERIOR Width x Depth x Height | TEMP RANGE (SETPOINT) | DOORS | SHELVES | DEFROST | CERTIFICATION | SHIPPING WT. kg / LBS. |
|-----------------|---------------------------------|---|--|-----------------------------|-------|---------|---------|---------------|---------------------------|
| K212GDR | 14/400 | W-480mm (18.9") D-500mm (19.69") H- 1150mm (45.28") | W-600mm (23.62") D-642mm (25.27") H-183mm (71.75") | 1 TO 12°C Preset +4°C | 1 | 4 | AUTO | UL / cUL | 75kg (165lbs) |