


K205ULT Ultra-low Temp Freezer

USES: / CLINICAL APPLICATIONS / VACCINES & MEDICATIONS / BIOLOGICS / DIAGNOSTICS KITS / MOLECULAR BIOLOGY / INDUSTRIAL TESTING

K2 Scientific's K205ULT -86°C freezer is designed for the most rigorous commercial applications. It has a single-compressor system built to run efficiently providing rapid cool down after door opening. The controls are easy to operate and feature temperature alarms.



- Intelligent microprocessor digital display, password protected
- Cabinet and ambient air temperature displays
- Adjustable temperature control
- Audible and visual high or low temperature alarms
- Power failure audible alarm
- Low-battery audible alarm
- Sensor failure alarm
- Key lock
- Interior cabinet constructed of stainless steel
- PUF-VIP, high density insulation
- Casters
- Remote, plug-in dry alarm contacts
- 0.5" temperature probe and data logger access port
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- Warranty: two-year full replacement for parts or compressor failure



MODEL NUMBER	CAPACITY CU. FT. / liters	EXTERIOR Height x Width x Depth	STORAGE AREA Width x Length x Depth	TEMP RANGE (SETPOINT)	DOORS	SHELVES	DEFROST	SHIPPING WT. LBS.
K205ULT	5 / 140	32.25" x 39.25" x 28.75"	18" x 25" x 18"	-50°C to -86°C	1	0	Manual	275 lbs