



HydraFlock® Flocked Swab

Ultrafine Tip with Polystyrene Handle - 100 mm breakpoint



HydraFlock®

Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patent Information: <https://www.puritanmedproducts.com/patents>

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.14" Dia. (3.556 mm)
- 0.66" Length (16.76 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 6" Length (152.4 mm)

Qualities

- CE marked, medical grade quality
- Ultrafine, highly absorbent flock tip
- Soft smooth tip
- Made in USA

Applications

- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection

Ref #	Packaging
3317-H Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3317-H Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)
25-3317-H 20 Sterile	50/20/1 (contains 1,000 individually wrapped swabs presented in 50 perforated strips of 20)