

SO-LOW CHEST STYLE FREEZER TO -40°C

MODEL C40-21

21 Cubic Ft. (594 liters) Chest – Temperature Range: 0°C to -40°C

Top Opening

Specifications:

| | |
|----------------------|---|
| * Capacity | 21 cubic ft. / 594 liters |
| * Temperature Range | 0°C to -40°C |
| * Dimensions | |
| * Exterior W x D x H | 91" x 36" x 49" 231 x 91 x 124 (cm) |
| * Interior W x D x H | 81" x 21.5" x 21" 205 x 54 x 53 (cm) |

TEMPERATURE CONTROL

Microprocessor control displays set point and chamber temperature.

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressor. CFC & HCFC free refrigerants.

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

115 volts 60 hertz, 1 phase.
20 AMP dedicated circuit required.
ETL electrical approval.

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.

