

SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-70RF

Refrigerator 46.6 cu. ft. -Freezer 24 cu. ft.

Upright – Temperature Range: 4°C (R) / -20°C (F) Three Doors

Specifications:

* Total Capacity		70 cubic ft.
* Temperature (Fridge)		4°C
* Temperature (Freezer)		-20°C
* Dimensions		
* Exterior w x d x h		78 x 34.75 x 81.5" 198 x 88 x 207 (cm)
* Fridge Interior w x d x h		48 x 28 x 60"
* Freezer Interior w x d x h		28 x 28 x 20"

Features:

Control	Digital Display with Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory For each compartment
Defrost	Auto
Doors	Three Solid, Locking
Casters	Four, Lockable
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Shelves	Epoxy Coated / 3 Per Compartment
Insulation:	Polyurethane
Evaporator	Fin and tube
Compressor	1/3 hp - 1/2 hp
Voltage:	115/60/1
Running Amps/Breaker	9.4/15 – 10.6/15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	950 lbs.
Warranty-(Domestic)	3 year parts and labor Compressor additional 4 years, no labor.



Specifications subject to change without notice