

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-15GD

15 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass Door

Specifications:

* Capacity	15 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	29.5" x 26.25" x 81" 74 x 66 x 205 (cm)
* Interior w x d x h	25.5" x 21.25" x 50.5" 64 x 53 x 125 (cm)

Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior finish	White powder coat
Casters	Four, locking
Shelves	Five
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	7.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice