

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-27SDSS

**27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Solid Door**

Specifications:

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	30" x 33.75" x 81" 76 x 85 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Door	Single solid, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, locking
Shelves	Three, Adjustable
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R290A CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0 /15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

