

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-30GD

30 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Two Glass Doors

Specifications:

* Capacity	30 cubic ft. / 849 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	47.25" x 30.25" x 81" 120 x 76 x 205 (cm)
* Interior w x d x h	43 x 22" x 54" 109 x 55 x 137 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Double glass – Lockable
Casters	Four, locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane
Refrigerant:	R290A
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps	8.3/15
Supply Plug	5-15 NEMA UL listed
Shipping Weight	500 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice