

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-50PT

50 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass doors on front and back

Specifications:

* Capacity	50 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	52" x 37.5" x 81.5" 132 x 95 x 207 (cm)
* Interior w x d x h	48" x 30" x 59" 121 x 76 x 149 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm Alarm Relay Dry Contacts Min/Max Memory
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Door	Glass doors front and back, lockable
Casters	Four 3" swivel, locking
Shelves	Six, Adjustable
Interior Light	One
Insulation:	Polyurethane foam
Refrigerant:	R-290A
Compressor	3/4 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage:	115/60/1
Running Amps	8.0
Supply Plug	5-15P NEMA UL listed
Shipping Weight	650 lbs.
Warranty-(Domestic)	3 year parts and labor. Compressors additional 4 years, no labor.



Specifications subject to change without notice.