

## *SO-LOW COUNTER TOP / FREE STANDING FREEZER MODEL MV20-2UCF*

**2.5 Cubic Ft. Upright – Temperature Range: 0°C to -25°C**

**One Solid Door**

### **Specifications:**

|                      |  |
|----------------------|--|
| * Capacity           | 2.5 cubic ft.                          |
| * Temperature Range  | 0°C to -25°C                           |
| * Dimensions         |  |
| * Exterior w x d x h | 20" x 21.5" x 25"<br>50 x 54 x 63 (cm) |
| * Interior w x d x h | 17" x 15" x 17"<br>43 x 38 x 43 (cm)   |

### **Features:**

|                       |  |
|-----------------------|--|
| Temperature Control   | Digital Display<br>With Hi/Low Alarm<br>Alarm Relay Dry Contacts |
| Exterior Construction | Painted White  |
| Access Port           | 1/2", included   |
| Door                  | Single Solid, locking  |
| Shelves               | Two full   |
| Insulation            | Polyurethane foam  |
| Refrigerant:          | R600   |
| Air Circulation       | Forced air   |
| Defrost               | Manual   |
| Voltage: Domestic     | 115/60/1   |
| Running Amps/Breaker  | 1.2/15   |
| Supply Plug           | 5-15P NEMA   |
| Shipping Weight       | 110 lbs.   |
| Warranty-(Domestic)   | Two year parts and labor   |



Specifications subject to change without notice