

## *SO-LOW COUNTER TOP / FREE STANDING REFRIGERATOR MODEL MV4-2UCRGDDA*

**2.5 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**One Glass Door**

### *Specifications:*

* Capacity	2.5 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	20" x 21.5" x 25" 50 x 54 x 63 (cm)
* Interior w x d x h	17" x 15" x 17" 43 x 38 x 43 (cm)

### *Features:*

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Exterior Construction	Painted White
Access Port	1/2"
Door	Single Glass, locking
Shelves	Two
Insulation	Polyurethane foam
Refrigerant:	R600
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	1.2/15
Supply Plug	5-15P NEMA
Shipping Weight	110 lbs.
Warranty-(Domestic)	Two year parts and labor



Specifications subject to change without notice