

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW UNDERCOUNTER REFRIGERATOR MODEL MV4-6UCRDA

4.0 Cubic Ft. Upright – Temperature Range: 2°C to 8°C

One Solid Door

Specifications:

* Capacity	4.0 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	23.75" x 21.5" x 32" 60 x 64 x 84 (cm)
* Interior w x d x h	20" x 15.25" x 22.75" 51 x 38 x 57 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2"
Door	single solid, locking
Shelves	Two adjustable, 1 full, 1 half
Insulation	Polyurethane foam
Refrigerant:	R600
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	1.2/15
Supply Plug	5-15P NEMA
Shipping Weight	150 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice