

SO-LOW UNDERCOUNTER FREEZER MODEL MV40-4UCF

4 Cubic Ft. Upright – Temperature Range: -30°C to -40°C

One Solid Door

Specifications:

| | |
|----------------------|---|
| * Capacity | 4 cubic ft. |
| * Temperature Range | -30°C to -40°C |
| * Dimensions | |
| * Exterior w x d x h | 24" x 24" x 33" 61 x 61 x 83 (cm) |
| * Interior w x d x h | 17.5" x 15.5" x 24.75" 44 x 39 x 62 (cm) |

Features:

| | |
|----------------------|--|
| Temperature Control | Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts |
| Access Port | 1/2" Rear |
| Door | Single solid, locking |
| Shelves | Three adjustable |
| Insulation | Polyurethane foam |
| Refrigerant: | R290 |
| Compressor | 1/4 hp |
| Defrost | Manual |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 5/15 |
| Supply Plug | 5-15P NEMA |
| Shipping Weight | 150 lbs. |
| Warranty-(Domestic) | Two year parts and labor |



Specifications subject to change without notice