

SO-LOW UNDERCOUNTER FREEZER MODEL MV85-1

**0.25 Cubic Ft. (7 liters) Upright – Temperature Range: -40°C to -85°C
 Solid Door**

Specifications:

* Capacity	0.25 cubic ft. / 7 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior w x d x h	15.75" x 25" x 26.25" 40 x 60 x 66.5 (cm)
* Interior w x d x h	5.75" x 5.5" x 12" 15 x 14.3 x 31 (cm)

Features:

Control	Digital Datalogging Battery Back-up Hi/Low Alarm
Clearance	4" on sides and back
Doors	Single Solid, locking
Legs	Four, leveling
Insulation	Polyurethane CFC Free
Evaporator	Cold wall
Defrost	Manual
Voltage: Domestic	115/60/1 – 2.8 amps
Supply Plug	NEMA 5-15P ETL listed
Shipping Weight	210 lbs.
Warranty-(Domestic)	One year parts and labor



Specifications subject to change without notice