

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## ***SO-LOW UNDERCOUNTER FREEZER MODEL MV85-2***

**2 Cubic Ft. (54 liters) Upright – Temperature Range: -40°C to -85°C  
Solid Door**

### ***Specifications:***

* Capacity	2 cubic ft. / 54 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior w x d x h	23" x 27" x 32" 60 x 70 x 81 (cm)
* Interior w x d x h	14" x 18" x 12.5" 36 x 46 x 32 (cm)

### ***Features:***

Control	Digital Datalogging Battery Back-up Hi/Low Alarm
Clearance	4" on sides and back
Doors	Single Solid, locking
Legs	Four, leveling
Insulation	Polyurethane CFC Free
Evaporator	Cold wall
Defrost	Manual
Voltage: Domestic	115/60/1 – 6.3 amps
Supply Plug	NEMA 5-15P ETL listed
Shipping Weight	210 lbs.
Warranty	One year parts and labor



Specifications subject to change without notice

### ***Accessories:***

Inventory Rack #17-15-2  
Rack holds Fifteen 2" high x 5.25" x 5.25"  
boxes with dividers. MV85-2 holds 2 racks.