

SO-LOW CHEST STYLE FREEZER TO -40°C

MODEL NC40-17

17 Cubic Ft. (480 liters) Chest – Temperature Range: 0°C to -40°C

Top Opening

Specifications:

* Capacity	17 cubic ft. / 480 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	81.25" x 35.5" x 48.5" 206 x 85 x 123 (cm)
* Interior W x D x H	69" x 20" x 21" 175 x 50 x 53 (cm)

COLOR TOUCH SCREEN TEMPERATURE CONTROL

Touch Screen control
Displays set point and chamber temperature
Hi/Low Alarm W/ Email and SMS notification
Alarm Relay Dry Contacts
Power Failure Dry Contacts
Battery Back-up
Viewable Temperature graph
Data logging downloadable via USB or FTP
Multi-Level Security
VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressors. CFC & HCFC free refrigerants

VOLTAGE

20 AMP dedicated circuit required.
ETL electrical approval.

Specifications subject to change without notice

