

SO-LOW CHEST STYLE FREEZER TO -40°C

MODEL NC40-5

5 Cubic Ft. (142 liters) Chest – Temperature Range: 0°C to -40°C
Top Opening

Specifications:

* Capacity	5 cubic ft. / 142 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	40.5" x 31" x 47" 102 x 81 x 119 (cm)
* Interior W x D x H	30" x 18" x 16" 76 x 46 x 40 (cm)

COLOR TOUCH SCREEN **TEMPERATURE CONTROL**

Touch Screen control
Displays set point and chamber temperature
Hi/Low Alarm W/ Email and SMS notification
Alarm Relay Dry Contacts
Power Failure Dry Contacts
Battery Back-up
Viewable Temperature graph
Data logging downloadable via USB or FTP
Multi-Level Security
VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressors. CFC & HCFC free refrigerants

VOLTAGE

15 AMP dedicated circuit required.
ETL electrical approval.

