

SO-LOW CHEST STYLE FREEZER TO -40°C

MODEL NC40-9

9 Cubic Ft. (255 liters) Chest – Temperature Range: 0°C to -40°C

Top Opening

Specifications:

| | |
|----------------------|--|
| * Capacity | 9 cubic ft. / 255 liters |
| * Temperature Range | 0°C to -40°C |
| * Dimensions | |
| * Exterior W x D x H | 57.5" x 29" x 46.5" 146 x 73 x 116 (cm) |
| * Interior W x D x H | 45.5" x 16" x 20" 115 x 40 x 50 (cm) |

COLOR TOUCH SCREEN TEMPERATURE CONTROL

Touch Screen control
Displays set point and chamber temperature
Hi/Low Alarm W/ Email and SMS notification
Alarm Relay Dry Contacts
Power Failure Dry Contacts
Battery Back-up
Viewable Temperature graph
Data logging downloadable via USB or FTP
Multi-Level Security
VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressors. CFC & HCFC free refrigerants

VOLTAGE

20 AMP dedicated circuit required.
ETL electrical approval.

Specifications subject to change without notice

