

SO-LOW CHEST STYLE FREEZER TO -85°C

MODEL NC85-12

12 Cubic Ft. (340 liters) Chest – Temperature Range: -40°C to -85°C

Top Opening

Specifications:

* Capacity	12 cubic ft. / 340 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior W x D x H	65.5" x 32" x 47.5" 166 x 81 x 120 (cm)
* Interior W x D x H	55" x 19" x 20" 139 x 48 x 50 (cm)

COLOR TOUCH SCREEN TEMPERATURE CONTROL

Touch Screen control
Displays set point and chamber temperature
Hi/Low Alarm W/ Email and SMS notification
Alarm Relay Dry Contacts
Power Failure Dry Contacts
Battery Back-up
Viewable Temperature graph
Data logging downloadable via USB or FTP
Multi-Level Security
VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

VOLTAGE

20 AMP dedicated circuit required.
ETL electrical approval.

Specifications subject to change without notice

