

SO-LOW CHEST STYLE FREEZER TO -85°C ***MODEL NC85-2***

2 Cubic Ft. (56 liters) Chest – Temperature Range: -40°C to -85°C
Mini Chest Freezer – Space Saver

Specifications:

| | |
|----------------------|---|
| * Capacity | 2 cubic ft. / 56 liters |
| * Temperature Range | -40°C to -85°C |
| * Dimensions | |
| * Exterior W x D x H | 25.25" x 25.75" x 40" 64 x 65 x 101 (cm) |
| * Interior W x D x H | 17" x 14" x 14" 43 x 35 x 35 (cm) |

COLOR TOUCH SCREEN **TEMPERATURE CONTROL**

Touch Screen control
Displays set point and chamber temperature
Hi/Low Alarm W/ Email and SMS notification
Alarm Relay Dry Contacts
Power Failure Dry Contacts
Battery Back-up
Viewable Temperature graph
Data logging downloadable via USB or FTP
Multi-Level Security
VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

VOLTAGE

115 volts 60 hertz, 1 phase.
15 AMP dedicated circuit required.
ETL electrical approval.

