

## *SO-LOW UPRIGHT FREEZER TO -40°C*

### *MODEL NU40-22*

**22 Cubic Ft. (626 liters) Upright – Temperature Range: 0°C to -40°C**

#### **Upright Freezer**

#### ***Specifications:***

* Capacity	22 cubic ft. / 626 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	43" x 37" x 79.5" 109 x 93 x 201 (cm)
* Interior W x D x H	28.5" x 26" x 51.25" 72 x 66 x 130 (cm)

#### **COLOR TOUCH SCREEN TEMPERATURE CONTROL**

Touch Screen control  
Displays set point and chamber temperature  
Hi/Low Alarm W/ Email and SMS notification  
Alarm Relay Dry Contacts  
Power Failure Dry Contacts  
Battery Back-up  
Viewable Temperature graph  
Data logging downloadable via USB or FTP  
Multi-Level Security  
VNC Remote Access

#### **ALARM SYSTEM WITH RELAY**

The battery operated alarm system will emit an audible & visual signal when there is a mechanical or electrical failure. Relay for remote alarm hook-up.

#### **REFRIGERATION SYSTEM**

A So-Low single stage system with one hermetic compressor. CFC & HCFC free refrigerants.

#### **CONSTRUCTION**

Chamber is 14-ga. Zinc coated galvanized steel.  
Exterior is 16-ga. steel. Power coated cool gray.

#### **VOLTAGE**

20 AMP dedicated circuit required.  
ETL electrical approval.

