

## *SO-LOW UPRIGHT FREEZER TO -40°C*

### *MODEL U40-25*

**25 Cubic Ft. (710 liters) Upright – Temperature Range: 0°C to -40°C**

#### **Upright Freezer**

#### ***Specifications:***

* Capacity	25 cubic ft. / 710 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	50" x 37" x 79.5" 127 x 93 x 201 (cm)
* Interior W x D x H	34" x 25" x 51.25" 86 x 63 x 130 (cm)

#### **TEMPERATURE CONTROL**

Microprocessor control displays set point and chamber temperature.

#### **ALARM SYSTEM WITH RELAY**

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

#### **REFRIGERATION SYSTEM**

A So-Low single stage system with one hermetic compressor. CFC & HCFC free refrigerants.

#### **CONSTRUCTION**

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

#### **VOLTAGE**

20 AMP dedicated circuit required.  
ETL electrical approval.

#### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

Specifications subject to change without notice

